



medw|ork®



POL_{Series}
Polypectomy Snares

POL_{Series} Polypectomy Snares




With its use of unique materials combined with intelligent design, the POL series set the standard in polypectomy:

The snare loop can be manipulated with a very high degree of precision, facilitating precise positioning and the complete resection of polyps, especially in difficult anatomical positions.

Polypectomy snares with integrated sclerotherapy needle (POL1-H3 . . .) enable you to safely elevate the tissue layers to be ablated, especially when carrying out mucosal resection.

POL^{Series} Polypectomy Snares



Compacted Wire Technology
Reliable coagulation with clean cut




Combination of flexibility and rigidity
Precise control of the snare loop

POL1 _ Disposable polypectomy snares _ Oval

Product code Order number	Shape	Characteristics	ø Snare	ø Tube	Length	Handle	Unit
POL1-B1-10-23-220-OL 501891	Oval	Braided wire, Compacted Wire Technology	10 mm	2.3 mm	220 cm	OL power connection	10
POL1-B1-15-23-220-OL 501207	Oval	Braided wire, Compacted Wire Technology	15 mm	2.3 mm	220 cm	OL power connection	10
POL1-B1-30-23-220-OL 501208	Oval	Braided wire, Compacted Wire Technology	30 mm	2.3 mm	220 cm	OL power connection	10
POL1-B1-40-23-220-OL 501210	Oval	Braided wire, Compacted Wire Technology	40 mm	2.3 mm	220 cm	OL power connection	10
POL1-B9-15-23-220-OL 502050	Oval	Monofil	15 mm	2.3 mm	220 cm	OL power connection	10
POL1-B9-22-23-220-OL 502048	Oval	Monofil	22 mm	2.3 mm	220 cm	OL power connection	10
POL1-B9-30-23-220-OL 502049	Oval	Monofil	30 mm	2.3 mm	220 cm	OL power connection	10





POL_{Series} Polypectomy Snares



POL1 _ Disposable polypectomy snares _ Oval, rotatable

Product code Order number	Shape	Characteristics	ø Snare	ø Tube	Length	Handle	Unit
POL1-B3-10-23-220-OL 502029	Oval	Monofil, rotatable	10 mm	2.3 mm	220 cm	OL power connection	10
POL1-B3-15-23-220-OL 501215	Oval	Monofil, rotatable	15 mm	2.3 mm	220 cm	OL power connection	10
POL1-B3-22-23-220-OL 501216	Oval	Monofil, rotatable	22 mm	2.3 mm	220 cm	OL power connection	10
POL1-B3-30-23-220-OL 501217	Oval	Monofil, rotatable	30 mm	2.3 mm	220 cm	OL power connection	10
POL1-B8-22-23-220-OL 501219	Oval	Braided wire, rotatable	22 mm	2.3 mm	220 cm	OL power connection	10
POL1-B8-30-23-220-OL 501220	Oval	Braided wire, rotatable	30 mm	2.3 mm	220 cm	OL power connection	10




POL^{Series} Polypectomy Snares



The right solution for every challenge


POL1 _ Disposable polypectomy snares _ Hexagonal

Product code Order number	Shape	Characteristics	ø Snare	ø Tube	Length	Handle	Unit
POL1-F1-30-23-220-OL 502101	Hexagonal	Braided wire	30 mm	2.3 mm	220 cm	OL power connection	10




POL1 _ Disposable polypectomy snares _ Asymmetrical, combi

Product code Order number	Shape	Characteristics	ø Snare	ø Tube	Length	Handle	Unit
POL1-K1-25-23-220-OL 502102	Asymmetrical, combi	Braided wire, Compacted Wire Technology	15 / 25 mm	2.3 mm	220 cm	OL power connection	10



POL1 _ Disposable polypectomy snares _ Oval, economy

Product code Order number	Shape	Characteristics	ø Snare	ø Tube	Length	Handle	Unit
POL1-B7-20-23-220-OL 501903	Oval	Braided wire	20 mm	2.3 mm	220 cm	OL power connection	10
POL1-B7-30-23-220-OL 501218	Oval	Braided wire	30 mm	2.3 mm	220 cm	OL power connection	10



POL^{Series} Polypectomy Snares




Resection Master _ Disposable polypectomy snares for mucosal resection

Product code Order number	Shape	Characteristics	ø Snare	ø Tube	Length	Handle	Unit
POL1-C5-22-23-220-OL 501896	Oval	Flat metal band, proximal strengthened	22 mm	2.3 mm	220 cm	OL power connection	10
POL1-C5-30-23-220-OL 501222	Oval	Flat metal band, proximal strengthened	30 mm	2.3 mm	220 cm	OL power connection	10




POL^{Series} Polypectomy Snares

Integrated needle
Optimum cushion formation



Snare Inflator _ Combi instrument

Product code Order number	Shape	Characteristics	ø Snare	ø Tube	Length	Handle	Unit
POL1-H3-22-27-230-OL 501231	Combi instrument	Polypectomy snare and injection needle in one instrument, monofil	22 mm	2.7 mm	230 cm	OL power connection	5
POL1-H3-30-27-230-OL 501948	Combi instrument	Polypectomy snare and injection needle in one instrument, monofil	30 mm	2.7 mm	230 cm	OL power connection	5



See for yourself the benefits of the

POL Series

Our national medical products consultants and international partners are happy to advise you in detail.

Contact us!



german||made