



# BAL

Series

**Stone Extraction Balloons**

# BAL<sup>Series</sup> Stone Extraction Balloons



## Stone Extraction Balloons

**BAL Series balloon catheters** are instruments noted for their very high quality and strength. The balloon strength also makes it ideal for use for stones with sharp edges or points.



### 2-lumen balloon catheter: BAL1-B2-16-75-200-35<sup>®</sup>









#### Benefit from the features:


- \_ accurate injection of contrast agent with the guide wire in place
- \_ tapered distal catheter from 7 to 5 French
- \_ excellent depth orientation due to markings on the catheter
- \_ x-ray markers both distally and proximally to the balloon
- \_ extremely strong balloon

# BAL Series Stone Extraction Balloons

## BAL1 \_ Disposable stone extraction balloons




Product Code <i>Order Number</i>	Shape	Characteristics	∅ Balloon	∅ Tube	Length	Accepted Guide Wire	Unit
<b>BAL1-C2-10-05-200-21</b> 500009	 Distal and proximal radiopaque markings	Distal contrast flow, 2 lumen	 10 mm	5 F	200 cm	.021"	2
<b>BAL1-B2-13-07-200-35</b> 500004	 Distal and proximal radiopaque markings	Distal contrast flow, 2 lumen	 13 mm	7 F	200 cm	.035"	2
<b>BAL1-B2-16-75-200-35</b> 500005	Distal and proximal radiopaque markings, tapered	Distal contrast flow, 2 lumen, <b>with sealing system</b>	 16 mm	7/5 F	200 cm	.035"	2
<b>BAL1-B2-16-76-200-35</b> 501926	Distal and proximal radiopaque markings, tapered	Distal contrast flow, 2 lumen	 16 mm	7/5 F	200 cm	.035"	10
<b>BAL1-C2-16-07-200-35</b> 500010	Distal and proximal radiopaque markings	Distal contrast flow, 2 lumen	 16 mm	7 F	200 cm	.035"	2
<b>BAL1-C3-18-07-200-35</b> 500011	Distal and proximal radiopaque markings	Distal contrast flow, 3 lumen	 18 mm	7 F	200 cm	.035"	2



## BAL1 \_ Disposable **easySure** stone extraction balloons



Product Code <i>Order Number</i>	Shape	Characteristics	∅ Balloon	∅ Tube	Length	Accepted Guide Wire	Unit
<b>BAL1-C3-18-07-200E35</b> 501916	Distal and proximal radiopaque markings	Distal contrast flow, 3 lumen	 18 mm	7 F	200 cm	.035"	2



 = Marked products may have a delivery time of up to 8 weeks.

Errors excepted.