




BAS

Series

Stone Extraction Baskets

BAS^{Series} Stone Extraction Baskets



Stone Extraction baskets

The perfect instrument for stone retrieval with the BAS Series stone extraction baskets. High-quality wires with outstanding expansion force provide excellent handling in the most difficult situations.



Double-lumen guide-wire stone extraction basket: BAS1-H1-20-27-200-35®

Particularly suitable for removal of intrahepatic stones and stones in the pancreatic duct


- _ easy access over a guide wire
- _ suitable for emergency lithotripsy
- _ excellent guide wire control
- _ access even with demanding stenoses

BAS Series Stone Extraction Baskets

BAS1 _ Disposable stone extraction baskets _ Spherical („twist n´CATCH“)



Product Code Order Number	Shape	Characteristics	ø Basket	ø Tube	Length	Unit
BAS1-K2-25-23-200 500027	Spherical	Braided wire, 0,4 mm, 130° rotation	25 mm	2,3 mm	200 cm	2



BAS1 _ Disposable guide wire stone extraction baskets _ 0,3 mm braided wire



Product Code Order Number	Shape	Characteristics	ø Basket	ø Tube	Length	Accepted Guide Wire	Unit
BAS1-H1-15-27-200-35 502018	Guide wire basket, half twisted	Braided wire, 0,3 mm	15 mm	2,7 mm	200 cm	.035"	2
BAS1-H1-20-27-200-35 502019	Guide wire basket, half twisted	Braided wire, 0,3 mm	20 mm	2,7 mm	200 cm	.035"	2
BAS1-H1-25-27-200-35 502020	Guide wire basket, half twisted	Braided wire, 0,3 mm	25 mm	2,7 mm	200 cm	.035"	2



BAS1 _ Disposable **easyCure** guide wire stone extraction baskets _ 0,3 mm braided wire




Product Code Order Number	Shape	Characteristics	ø Basket	ø Tube	Length	Accepted Guide Wire	Unit
BAS1-H1-20-27-200E35 501913	Guide wire basket, half twisted	Braided wire, 0,3 mm	20 mm	2,7 mm	200 cm	.035"	2



BAS1 _ Disposable stone extraction baskets _ 0,35 mm braided wire



Product Code Order Number	Shape	Characteristics	ø Basket	ø Tube	Length	Unit
BAS1-E2-20-23-200 500016	Fully twisted, 6 braids	Braided wire, 0,35 mm	20 mm	2,3 mm	200 cm	2



medwork® does not accept liability for system components of the BAS Series that are used with system components from other manufacturers.
For detailed information please see the instruction manual.

Errors excepted.

BAS Series Stone Extraction Baskets

BAS1 _ Disposable stone extraction baskets _ 0,4 mm braided wire



Product Code <i>Order Number</i>	Shape	Characteristics	ø Basket	ø Tube	Length	Unit
BAS1-A2-20-23-200 500018	Hexagonal	Braided wire, 0,4 mm	20 mm	2,3 mm	200 cm	2
BAS1-A2-25-23-200 500020	Hexagonal	Braided wire, 0,4 mm	25 mm	2,3 mm	200 cm	2
BAS1-A2-30-23-200 500022	Hexagonal	Braided wire, 0,4 mm	30 mm	2,3 mm	200 cm	2
BAS1-A2-35-23-200 500024	Hexagonal	Braided wire, 0,4 mm	35 mm	2,3 mm	200 cm	2



BAS1 _ Disposable stone extraction baskets for bronchoscopy



Product Code <i>Order Number</i>	Shape	Characteristics	ø Basket	ø Tube	Length	Unit
BAS1-C1-20-18-120 502010	Half twisted	Braided wire, 0,3 mm	20 mm	1,8 mm	120 cm	2



BAS1 _ Disposable stone extraction baskets for balloon enteroscopy



Product Code <i>Order Number</i>	Shape	Characteristics	ø Basket	ø Tube	Length	Unit
BAS1-C1-20-18-260 500025	Half twisted	Braided wire, 0,3 mm	20 mm	1,8 mm	260 cm	2
BAS1-C2-25-23-260 500026	Half twisted	Braided wire, 0,4 mm	25 mm	2,3 mm	260 cm	2




medwork® does not accept liability for system components of the BAS Series that are used with system components from other manufacturers. For detailed information please see the instruction manual.

Errors excepted.

BAS^{Series} Stone Extraction Baskets




BAS2 _ Reusable stone extraction baskets _ Spherical („twist'n´CATCH“)

134°C

Product Code <i>Order Number</i>	Shape	Characteristics	ø Basket	ø Tube	Length	Unit	
BAS2-K2-25-23-200 500120	Spherical	Braided wire, 0,4 mm, 130° rotation	25 mm	2,3 mm	200 cm	1	


BAS2 _ Reusable stone extraction baskets _ 0,3 mm braided wire

134°C

Product Code <i>Order Number</i>	Shape	Characteristics	ø Basket	ø Tube	Length	Unit	
BAS2-C1-15-18-200 500096	 Half twisted	Braided wire, 0,3 mm	15 mm	1,8 mm	200 cm	1	
BAS2-C1-20-18-200 500098	 Half twisted	Braided wire, 0,3 mm	20 mm	1,8 mm	200 cm	1	


BAS2 _ Reusable stone extraction baskets _ 0,35 mm braided wire

134°C


Product Code <i>Order Number</i>	Shape	Characteristics	ø Basket	ø Tube	Length	Unit	
BAS2-E2-30-23-200 500110	Fully twisted, 6 braids	Braided wire, 0,35 mm	30 mm	2,3 mm	200 cm	1	

BAS2 _ Reusable stone extraction baskets _ 0,4 mm braided wire

134°C

Product Code <i>Order Number</i>	Shape	Characteristics	ø Basket	ø Tube	Length	Unit	
BAS2-A2-20-23-200 500076	Hexagonal	Braided wire, 0,4 mm	20 mm	2,3 mm	200 cm	1	
BAS2-A2-25-23-200 500083	Hexagonal	Braided wire, 0,4 mm	25 mm	2,3 mm	200 cm	1	
BAS2-A2-30-23-200 500084	Hexagonal	Braided wire, 0,4 mm	30 mm	2,3 mm	200 cm	1	
BAS2-A2-35-23-200 500085	Hexagonal	Braided wire, 0,4 mm	35 mm	2,3 mm	200 cm	1	

medwork® does not accept liability for system components of the BAS Series that are used with system components from other manufacturers.
For detailed information please see the instruction manual.


 = Marked products may have a delivery time of up to 8 weeks.

Errors excepted.

BAS Series Stone Extraction Baskets



BAS2 _ Reusable stone extraction baskets _ 0,45 mm braided wire

134°C

Product Code <i>Order Number</i>	Shape	Characteristics	ø Basket	ø Tube	Length	Unit	
BAS2-A3-20-23-200 500088	 Hexagonal	Braided wire, 0,45 mm	20 mm	2,3 mm	200 cm	1	
BAS2-A3-25-23-200 500089	 Hexagonal	Braided wire, 0,45 mm	25 mm	2,3 mm	200 cm	1	
BAS2-B3-30-23-200 500094	 Combi	Braided wire, 0,45 mm	15 / 30 mm	2,3 mm	200 cm	1	
BAS2-C3-20-23-200 500105	Half twisted	Braided wire, 0,45 mm	20 mm	2,3 mm	200 cm	1	
BAS2-C3-25-23-200 500106	Half twisted	Braided wire, 0,45 mm	25 mm	2,3 mm	200 cm	1	
BAS2-C3-30-23-200 500107	Half twisted	Braided wire, 0,45 mm	30 mm	2,3 mm	200 cm	1	
BAS2-C3-35-23-200 500108	Half twisted	Braided wire, 0,45 mm	35 mm	2,3 mm	200 cm	1	



BAS2 _ Reusable stone extraction baskets _ Nitinol (not suitable for lithotripsy)

134°C


Product Code <i>Order Number</i>	Shape	Characteristics	ø Basket	ø Tube	Length	Unit	
BAS2-A4-20-23-200 500093	Hexagonal	Nitinol	20 mm	2,3 mm	200 cm	1	
BAS2-F4-15-23-200 500112	Fully twisted, 8 braids	Nitinol	15 mm	2,3 mm	200 cm	1	

BAS2 _ Reusable stone extraction baskets for bronchoscopy

134°C

Product Code <i>Order Number</i>	Shape	Characteristics	ø Basket	ø Tube	Length	Unit	
BAS2-A1-20-18-120 500074	 Hexagonal	Braided wire, 0,3 mm	20 mm	1,8 mm	120 cm	1	

medwork® does not accept liability for system components of the BAS Series that are used with system components from other manufacturers.
For detailed information please see the instruction manual.

 = Marked products may have a delivery time of up to 8 weeks.

Errors excepted.