

# CAN

Series

**ERCP-Cannulas**

medw**ork**.®

CAN



## Atraumatic and safe

We have great experience in of tip shaping of ERCP cannulas: we offer the right tip for every situation.



## 2-lumen ERCP-cannula: CAN1-B4-18-200-35<sup>®</sup>

### Benefit from the features:

- \_ excellent guide wire control
- \_ the second lumen removes the problem of sticking a guidewire
- \_ filiform, atraumatic tip

# CAN Series ERCP-Cannulas

## CAN1 \_ Disposable ERCP-cannulas accepting .035" guide wires



Product Code <i>Order Number</i>	Shape	Characteristics	ø Tube	Length	Accepted Guide wire	Unit	
<b>CAN1-A3-18-200-35</b> 500337	Standard	Radiopaque marking, 1 lumen	1,8 mm	200 cm	.035"	10	
<b>CAN1-A5-18-200-35</b> 500339	Standard	Metal tip, 1 lumen	1,8 mm	200 cm	.035"	10	
<b>CAN1-B3-18-200-35</b> 500343	Tapered	Radiopaque marking, 1 lumen	1,8 mm	200 cm	.035"	10	
<b>CAN1-B4-18-200-35</b> 502008	Tapered	Radiopaque marking, <b>2 lumen</b>	1,8 mm	200 cm	.035"	10	
<b>CAN1-D5-18-200-35</b> 500348	Bottle neck	Metal tip, 1 lumen	1,8 mm	200 cm	.035"	10	

## CAN1 \_ Disposable ERCP-cannulas for balloon enteroscopy



Product Code <i>Order Number</i>	Shape	Characteristics	ø Tube	Length	Accepted Guide wire	Unit	
<b>CAN1-B3-18-260-35</b> 501856	Tapered	Radiopaque marking, 1 lumen	1,8 mm	260 cm	.035"	2	




## CAN1 \_ Disposable **easycure** ERCP-cannulas



Product Code <i>Order Number</i>	Shape	Characteristics	ø Tube	Length	Accepted Guide wire	Unit	
<b>CAN1-B4-18-200-E35</b> 502007	Tapered	Radiopaque marking, <b>2 lumen</b>	1,8 mm	200 cm	.035"	2	


## CAN2 \_ Reusable ERCP-cannulas accepting .018" guide wires

134°C

Product Code <i>Order Number</i>	Shape	Characteristics	ø Tube	Length	Accepted Guide wire	Unit	
<b>CAN2-E3-16-200-18</b> 500401	Long tapered tip	Radiopaque marking, 1 lumen	1,6 mm	200 cm	.018"	1	
<b>CAN2-F5-16-200-18</b> 500404	 Minor papilla	Metal tip, 1 lumen	1,6 mm	200 cm	.018"	1	






## CAN2 \_ Reusable ERCP-cannulas accepting .021" guide wires

134°C

Product Code <i>Order Number</i>	Shape	Characteristics	ø Tube	Length	Accepted Guide wire	Unit	
<b>CAN2-L5-16-200-21</b> 500407	Long, bottle neck	Metal tip, 1 lumen	1,6 mm	200 cm	.021"	1	

## CAN2 \_ Reusable ERCP-cannulas accepting .025" guide wires

134°C

Product Code <i>Order Number</i>	Shape	Characteristics	ø Tube	Length	Accepted Guide wire	Unit	
<b>CAN2-B1-16-200-25</b> 500377	Tapered	Standard, 1 lumen	1,6 mm	200 cm	.025"	1	
<b>CAN2-B5-16-200-25</b> 500384	 Tapered	Metal tip, 1 lumen	1,6 mm	200 cm	.025"	1	
<b>CAN2-C5-16-200-25</b> 500387	Conical	Metal tip, 1 lumen	1,6 mm	200 cm	.025"	1	
<b>CAN2-D5-16-200-25</b> 500396	Bottle neck	Metal tip, 1 lumen	1,6 mm	200 cm	.025"	1	










 = Marked products may have a delivery time of up to 8 weeks.

Errors excepted.

# CAN Series ERCP-Cannulas

CAN2 \_ Reusable ERCP-cannulas accepting .035" guide wires

134°C

Product Code <i>Order Number</i>	Shape	Characteristics	ø Tube	Length	Accepted Guide wire	Unit	
CAN2-A1-18-200-35 500368	Standard	Standard, 1 lumen	1,8 mm	200 cm	.035"	1	
CAN2-A2-18-200-35 500369	 Standard	Standard, Y-adapter, 1 lumen	1,8 mm	200 cm	.035"	1	
CAN2-A3-18-200-35 500371	 Standard	Radiopaque marking, 1 lumen	1,8 mm	200 cm	.035"	1	
CAN2-A4-18-200-35 500372	Standard	Radiopaque marking, Y-adapter, 1 lumen	1,8 mm	200 cm	.035"	1	
CAN2-A5-18-200-35 500374	Standard	Metal tip, 1 lumen	1,8 mm	200 cm	.035"	1	
CAN2-A6-18-200-35 500375	 Standard	Metal tip, Y-adapter, 1 lumen	1,8 mm	200 cm	.035"	1	
CAN2-B1-18-200-35 500378	Tapered	Standard, 1 lumen	1,8 mm	200 cm	.035"	1	
CAN2-B2-18-200-35 500379	 Tapered	Standard, Y-adapter, 1 lumen	1,8 mm	200 cm	.035"	1	
CAN2-B3-18-200-35 500381	Tapered	Radiopaque marking, 1 lumen	1,8 mm	200 cm	.035"	1	
CAN2-B4-18-200-35 500383	Tapered	Radiopaque marking, Y-adapter, 1 lumen	1,8 mm	200 cm	.035"	1	
CAN2-B5-18-200-35 500385	Tapered	Metal tip, 1 lumen	1,8 mm	200 cm	.035"	1	
CAN2-B6-18-200-35 500386	Tapered	Metal tip, Y-adapter, 1 lumen	1,8 mm	200 cm	.035"	1	
CAN2-C5-18-200-35 500388	Conical	Metal tip, 1 lumen	1,8 mm	200 cm	.035"	1	
CAN2-C6-18-200-35 500389	Conical	Metal tip, Y-adapter, 1 lumen	1,8 mm	200 cm	.035"	1	
CAN2-D5-18-200-35 500398	Bottle neck	Metal tip, 1 lumen	1,8 mm	200 cm	.035"	1	
CAN2-D6-18-200-35 500400	Bottle neck	Metal tip, Y-adapter, 1 lumen	1,8 mm	200 cm	.035"	1	

 = Marked products may have a delivery time of up to 8 weeks.

Errors excepted.