

# CLE

Series

**Cleaning Brushes**

medw**ork**.<sup>®</sup>

CLE

# CLE<sup>Series</sup> Cleaning Brushes



## Efficient and thorough

The CLE Series brushes provide what you are looking for: best possible hygienic safety and protection from bacterial transmission. Designs with two and three brush heads not only clean particularly effectively but also save time

### Combination cleaning brush: CLE1-C2-50-18-230<sup>®</sup>

#### Benefit from the features:

- \_ cleans all areas of the scope
- \_ low friction shaft
- \_ tapered reinforced tip

# CLE Series Cleaning Brushes

## CLE1 \_ Disposable cleaning brushes



Product Code <i>Order Number</i>	Shape	Characteristics	ø Bristle	ø Tube	Length	Unit
<b>CLE1-B1-30-16-100</b> <i>500458</i>	Bronchoscope	1 Bristle	3,0 mm	1,6 mm	100 cm	50
<b>CLE1-B2-30-16-100</b> <i>500459</i>	Bronchoscope	2 Bristles	3,0 mm	1,6 mm	100 cm	50
<b>CLE1-C1-50-18-230</b> <i>500460</i>	Colonoscope	1 Bristle	5,0 mm	1,8 mm	230 cm	50
<b>CLE1-C1-65-18-230</b> <i>500462</i>	Colonoscope	1 Bristle	6,5 mm	1,8 mm	230 cm	50
<b>CLE1-C2-50-18-230</b> <i>500464</i>	Colonoscope	2 Bristles	5,0 mm	1,8 mm	230 cm	50
<b>CLE1-C3-50-18-230</b> <i>500466</i>	Colonoscope	3 Bristles	5,0 mm	1,8 mm	230 cm	50
<b>CLE1-G1-30-16-180</b> <i>500468</i>	Gastroscope	1 Bristle	3,0 mm	1,6 mm	180 cm	50
<b>CLE1-G1-50-18-180</b> <i>500469</i>	Gastroscope	1 Bristle	5,0 mm	1,8 mm	180 cm	50
<b>CLE1-G2-50-18-180</b> <i>500471</i>	Gastroscope	2 Bristles	5,0 mm	1,8 mm	180 cm	50
<b>CLE1-V2-016</b> <i>500480</i>	Endoscope valve intake	2 Bristles	11 mm / 5 mm	-	16 cm	10



## CLE1 \_ Disposable cleaning brushes for balloon enteroscopy










Product Code <i>Order Number</i>	Shape	Characteristics	ø Bristle	ø Tube	Length	Unit
<b>CLE1-C2-30-16-280</b> <i>500463</i>	Colonoscope	2 Bristles	3,0 mm	1,6 mm	280 cm	10
<b>CLE1-C2-50-18-280</b> <i>500465</i>	Colonoscope	2 Bristles	5,0 mm	1,8 mm	280 cm	10



## CLE2 \_ Reusable cleaning brushes

134°C

Product Code <i>Order Number</i>	Shape	Characteristics	ø Bristle	ø Tube	Length	Unit	
<b>CLE2-C1-65-20-220</b> 500488	 Colonoscope	1 Bristle, distal protection ball	6,5 mm	2,0 mm	220 cm	1	
<b>CLE2-C2-65-20-220</b> 500491	 Colonoscope	2 Bristles, distal protection ball	6,5 mm	2,0 mm	220 cm	1	
<b>CLE2-G1-65-20-180</b> 500499	 Gastroscope	1 Bristle, distal protection ball	6,5 mm	2,0 mm	180 cm	1	
<b>CLE2-G2-65-20-180</b> 500501	 Gastroscope	2 Bristles, distal protection ball	6,5 mm	2,0 mm	180 cm	1	
<b>CLE2-O1-36-08-090</b> 500503	 Oesophageal dilator	1 Bristle, distal protection ball	3,6 mm	0,8 mm	90 cm	1	

 = Marked products may have a delivery time of up to 8 weeks.

Errors excepted.