


medw^ork[®]

DIL

DIL

Series

Dilation Devices



Flexible and atraumatic

The DIL series dilatation catheters and bougies atraumatically expand benign and malignant stenoses in the gastrointestinal tract. The very good radiological visibility makes safe application and control very easy.

3-Stage balloon dilators: DIL1-A1-06-060-18-35[®]


Benefit from the features:

- _ 3 sizes of dilation for each balloon
- _ pre-lubricated balloon shaft
- _ glow in the dark size and inflation information

DIL Series Dilation Devices


DIL1 _ Disposable dilation balloons for the gastrointestinal and bile tract (OTW / TTS) ②

Product Code <i>Order Number</i>	Shape	Characteristics	∅ Balloon	Length Balloon	∅ Tube	Accepted Guide Wire	Unit
DIL1-A1-06-060-18-35 500520	Dilation balloon	2 Radiopaque markings, length 180 cm	6 mm	60 mm	1,8 mm	.035"	1
DIL1-A1-08-060-18-35 500523	Dilation balloon	2 Radiopaque markings, length 180 cm	8 mm	60 mm	1,8 mm	.035"	1
DIL1-A1-12-040-23-35 500524	Dilation balloon	2 Radiopaque markings, length 180 cm	12 mm	40 mm	2,3 mm	.035"	1
DIL1-A1-15-040-23-35 500525	Dilation balloon	2 Radiopaque markings, length 210 cm	15 mm	40 mm	2,3 mm	.035"	1
DIL1-A1-18-040-23-35 500526	Dilation balloon	2 Radiopaque markings, length 210 cm	18 mm	40 mm	2,3 mm	.035"	1
DIL1-A1-20-045-23-35 500528	Dilation balloon	2 Radiopaque markings, length 210 cm	20 mm	45 mm	2,3 mm	.035"	1



















DIL4 _ Disposable insufflation manometer ②


Product Code <i>Order Number</i>	Characteristics	Unit
DIL4-A1-30 500561	With lockable piston and rotatable handle, suitable for: DIL1-A1-06-060-18-35 bis -DIL1-A1-12-040-23-35 measurement range up to 30 bar vacuum, scaling in bar and PSI	1



DIL2 _ Reusable oesophagus bougies (OTW)

Product Code <i>Order Number</i>	Shape	Characteristics	ø Bougie	Length	Accepted Guide Wire	Unit
DIL2-C2-05-800-35 500531	 Bougie	Completely radiopaque, outer markings	5 mm	800 mm	.035"	1
DIL2-C2-06-800-35 500533	 Bougie	Completely radiopaque, outer markings	6 mm	800 mm	.035"	1
DIL2-C2-07-800-35 500535	 Bougie	Completely radiopaque, outer markings	7 mm	800 mm	.035"	1
DIL2-C2-08-800-35 500537	 Bougie	Completely radiopaque, outer markings	8 mm	800 mm	.035"	1
DIL2-C2-09-800-35 500539	 Bougie	Completely radiopaque, outer markings	9 mm	800 mm	.035"	1
DIL2-C2-10-800-35 500541	 Bougie	Completely radiopaque, outer markings	10 mm	800 mm	.035"	1
DIL2-C2-11-800-35 500543	 Bougie	Completely radiopaque, outer markings	11 mm	800 mm	.035"	1
DIL2-C2-12-800-35 500545	 Bougie	Completely radiopaque, outer markings	12 mm	800 mm	.035"	1
DIL2-C2-13-800-35 500547	 Bougie	Completely radiopaque, outer markings	13 mm	800 mm	.035"	1
DIL2-C2-14-800-35 500548	 Bougie	Completely radiopaque, outer markings	14 mm	800 mm	.035"	1
DIL2-C2-15-800-35 500550	 Bougie	Completely radiopaque, outer markings	15 mm	800 mm	.035"	1
DIL2-C2-16-800-35 500552	 Bougie	Completely radiopaque, outer markings	16 mm	800 mm	.035"	1
DIL2-C2-17-800-35 500553	 Bougie	Completely radiopaque, outer markings	17 mm	800 mm	.035"	1
DIL2-C2-18-800-35 500555	 Bougie	Completely radiopaque, outer markings	18 mm	800 mm	.035"	1
DIL2-C2-19-800-35 500556	 Bougie	Completely radiopaque, outer markings	19 mm	800 mm	.035"	1
DIL2-C2-20-800-35 500558	 Bougie	Completely radiopaque, outer markings	20 mm	800 mm	.035"	1






 = Marked products may have a delivery time of up to 8 weeks.

Errors excepted.

DIL Series Dilation Devices

DIL3 _ Reusable oesophagus bougie sets (OTW)

Product Code <i>Order Number</i>	Characteristics	Unit
DIL3-C2-07-800 500559	 Consisting of: 7 Bougies (in 5, 7, 9, 11, 13, 14 and 15 mm), 1 guide wire, 1 cleaning brush, 1 carrying case	1
DIL3-C2-16-800 500560	 Consisting of: 16 Bougies (in 5 - 20 mm), 1 guide wire, 1 cleaning brush, 1 carrying brush	1

 = Marked products may have a delivery time of up to 8 weeks.

Errors excepted.