

medw<sup>o</sup>rk<sup>®</sup>

# INJ

Series  
**Injection Needles**

INJ



## Injection Needles

The high-quality lancet cut of the needle makes accurate injection easy with all products of the INJ Series.








## INJ series

### Benefit from the features:


- \_ the handle locking mechanism ensures it easy to feel that the needle is securely fixed in either the retracted or extended positions
- \_ protect the working channel of the endoscope with the distal metal hub with all INJ needles
- \_ all INJ Series injection needles are suitable for Histoacryl application

# INJ Series Injection Needles


## INJ1 \_ Disposable injection needles (click mechanism) ②

Product Code <i>Order Number</i>	Shape	Characteristics	ø Needle	Needle-lift	ø Tube	Length	Unit	
INJ1-A1-05-5-18-120 500744	Sclerosing	PTFE catheter, LL-connector, 1 lumen, atraumatic metal tip	0,5 mm	5 mm	1,8 mm	120 cm	10	
INJ1-A1-05-5-18-180 500745	Sclerosing	PTFE catheter, LL-connector, 1 lumen, atraumatic metal tip	0,5 mm	5 mm	1,8 mm	180 cm	10	
INJ1-A1-07-4-23-230 500750	Sclerosing	PTFE catheter, LL-connector, 1 lumen, atraumatic metal tip	0,7 mm	4 mm	2,3 mm	230 cm	10	
INJ1-A1-07-5-18-120 500751	Sclerosing	PTFE catheter, LL-connector, 1 lumen, atraumatic metal tip	0,7 mm	5 mm	1,8 mm	120 cm	10	
INJ1-A1-07-5-18-180 500752	Sclerosing	PTFE catheter, LL-connector, 1 lumen, atraumatic metal tip	0,7 mm	5 mm	1,8 mm	180 cm	10	
INJ1-A1-07-5-23-180 500757	Sclerosing	PTFE catheter, LL-connector, 1 lumen, atraumatic metal tip	0,7 mm	5 mm	2,3 mm	180 cm	10	
INJ1-A1-07-5-23-230 500760	Sclerosing	PTFE catheter, LL-connector, 1 lumen, atraumatic metal tip	0,7 mm	5 mm	2,3 mm	230 cm	10	
INJ1-A1-07-8-23-230 500771	Sclerosing	PTFE catheter, LL-connector, 1 lumen, atraumatic metal tip	0,7 mm	8 mm	2,3 mm	230 cm	10	
INJ1-A1-09-5-23-180 500772	Sclerosing	PTFE catheter, LL-connector, 1 lumen, atraumatic metal tip	0,9 mm	5 mm	2,3 mm	180 cm	10	
INJ1-A1-09-5-23-230 500773	Sclerosing	PTFE catheter, LL-connector, 1 lumen, atraumatic metal tip	0,9 mm	5 mm	2,3 mm	230 cm	10	

## INJ1 \_ Disposable injection needles for balloon enteroscopy (click mechanism) ②

Product Code <i>Order Number</i>	Shape	Characteristics	ø Needle	Needle-lift	ø Tube	Length	Unit	
INJ1-A1-07-5-18-260 500753	Sclerosing	PTFE catheter, LL-connector, 1 lumen, atraumatic metal tip	0,7 mm	5 mm	1,8 mm	260 cm	10	
INJ1-A1-07-5-23-260 500763	Sclerosing	PTFE catheter, LL-connector, 1 lumen, atraumatic metal tip	0,7 mm	5 mm	2,3 mm	260 cm	10	

## INJ1 \_ Disposable Fibrin needle (chamber technique) ②

Product Code <i>Order Number</i>	Shape	Characteristics	ø Needle	Needle-lift	ø Tube	Length	Unit	
INJ1-D1-08-5-26-210 500783	Fibrin / Thrombin	PTFE catheter, 2 LL-connectors, 2 lumen	0,8 mm	5 mm	2,6 mm	210 cm	2	

Errors excepted.