


medw<sup>o</sup>rk®

LIT

Series

Lithotripter Baskets

# LIT<sup>Series</sup> Lithotripter Baskets



## Lithotripter Baskets

**Ideal for any challenge that you face:** all products of this series are designed to allow you to replace the LIT basket with a guide wire at any time to cover particularly difficult stenoses.



### twist'n'CATCH: LIT1-K2-25-23-200<sup>®</sup>

#### Benefit from the features:


- \_ high expansion force with the special surface structure of the wire
- \_ 130° rotary movement while retracting the basket
- \_ infinite adjustment of different basket sizes
- \_ defined breaking point
- \_ atraumatic, filiform tip for simple passage through stenosis

# LIT Series Lithotripter Baskets

## LIT1 \_ Disposable spherical lithotripter basket \_ 0,4 mm braided wire („twist'n'CATCH“)




Product Code Order Number	Shape	Characteristics	ø Basket	ø Tube	Length	Unit
LIT1-K2-25-23-200 500826	Spherical	0,4 mm Braided wire, 130° rotation	25 mm	2,3 mm	200 cm / 400 cm	2



## LIT1 \_ Disposable spherical lithotripter basket for balloon enteroscopy \_ 0,4 mm braided wire („twist'n'CATCH“)




Product Code Order Number	Shape	Characteristics	ø Basket	ø Tube	Length	Unit
LIT1-K2-25-23-260 501904	Spherical	0,4 mm Braided wire, 130° rotation	25 mm	2,3 mm	260 cm / 400 cm	1



## LIT1 \_ Disposable lithotripter basket \_ 0,45 mm braided wire



Product Code Order Number	Shape	Characteristics	ø Basket	ø Tube	Length	Unit
LIT1-C3-20-23-200 500822	Half twisted	0,45 mm Braided wire	20 mm	2,3 mm	200 cm / 400 cm	2
LIT1-C3-25-23-200 500823	Half twisted	0,45 mm Braided wire	25 mm	2,3 mm	200 cm / 400 cm	2
LIT1-C3-30-23-200 500824	Half twisted	0,45 mm Braided wire	30 mm	2,3 mm	200 cm / 400 cm	2
LIT1-C3-35-23-200 500825	Half twisted	0,45 mm Braided wire	35 mm	2,3 mm	200 cm / 400 cm	2



## LIT2 \_ Reusable lithotripter basket \_ 0,3 mm braided wire

134°C

Product Code Order Number	Shape	Characteristics	ø Basket	ø Tube	Length	Unit
LIT2-C1-15-18-200 500871	Half twisted	0,3 mm Braided wire	15 mm	1,8 mm	200 cm / 400 cm	1
LIT2-C1-20-18-200 500872	Half twisted	0,3 mm Braided wire	20 mm	1,8 mm	200 cm / 400 cm	1




medwork® does not accept liability for system components of the LIT Series that are used with system components from other manufacturers. For detailed information please see the instruction manual.

Errors excepted.

## LIT2 \_ Reusable spherical lithotripter basket \_ 0,4 mm braided wire („twist'n'CATCH“)

134°C

Product Code Order Number	Shape	Characteristics	ø Basket	ø Tube	Length	Unit
LIT2-K2-25-23-200 500884	Spherical	0,4 mm Braided wire, 130° rotation	25 mm	2,3 mm	200 cm / 400 cm	1



## LIT2 \_ Reusable lithotripter basket \_ 0,45 mm braided wire

134°C

Product Code Order Number	Shape	Characteristics	ø Basket	ø Tube	Length	Unit
LIT2-A3-20-23-200 500866	Hexagonal	0,45 mm Braided wire	20 mm	2,3 mm	200 cm / 400 cm	1
LIT2-A3-25-23-200 500867	Hexagonal	0,45 mm Braided wire	25 mm	2,3 mm	200 cm / 400 cm	1
LIT2-A3-30-23-200 500868	Hexagonal	0,45 mm Braided wire	30 mm	2,3 mm	200 cm / 400 cm	1
LIT2-A3-35-23-200 500869	Hexagonal	0,45 mm Braided wire	35 mm	2,3 mm	200 cm / 400 cm	1
LIT2-B3-30-23-200 500870	Combi	0,45 mm Braided wire	15 / 30 mm	2,3 mm	200 cm / 400 cm	1
LIT2-C3-20-23-200 500877	Half twisted	0,45 mm Braided wire	20 mm	2,3 mm	200 cm / 400 cm	1
LIT2-C3-25-23-200 500878	Half twisted	0,45 mm Braided wire	25 mm	2,3 mm	200 cm / 400 cm	1
LIT2-C3-30-23-200 500879	Half twisted	0,45 mm Braided wire	30 mm	2,3 mm	200 cm / 400 cm	1
LIT2-C3-35-23-200 500880	Half twisted	0,45 mm Braided wire	35 mm	2,3 mm	200 cm / 400 cm	1







## LIT2 \_ Reusable lithotripter handle and metal springs

134°C

Product Code Order Number	Shape	Characteristics	Diameter	Length	Unit
LIT2-D-HAN 500881	-	Pull through technique	-	-	1
LIT2-G4-26-185 500882	Bile duct, transendoscopic	Spiral spring	2,6 mm	185 cm	1
LIT2-N4-38-090 500885	Bile duct, emergency	Spiral spring	3,8 mm	90 cm	1
LIT2-P4-23-185 500886	Pancreatic duct, transendoscopic	Spiral spring	2,3 mm	185 cm	1

medwork® does not accept liability for system components of the LIT Series that are used with system components from other manufacturers. For detailed information please see the instruction manual.

Errors excepted.