


medwörk®

PAP

Series

Papillotomes

PAP



Papillotomes

Expertise combined with technology: accurate control and perfect positioning, available in a wide range of models to offer you the best instrument for every situation.



3-lumen papillotome: PAP1-C7-30-35-OL[®]

Benefit from the features:

- _ insulated proximal cutting wire
- _ wire easily moved in while administering contrast agent simultaneously
- _ filoform curved distal end

PAP Series Papillotomes

PAP1 _ Disposable papillotomes _ 1 Lumen



Product Code <i>Order Number</i>	Shape	Characteristics	Cutting Length	Accepted Guide Wire	Handle	Unit	
PAP1-A1-20-OL 501027	Precut	1 Lumen, Braided wire	20 mm	-	OL power connection	2	
PAP1-C1-30-OL 501033	Tapered	1 Lumen, Braided wire	30 mm	-	OL power connection	2	
PAP1-K1-06-OL 501058	Needle	1 Lumen	6 mm	-	OL power connection	2	

PAP1 _ Disposable papillotomes _ 2 Lumen



Product Code <i>Order Number</i>	Shape	Characteristics	Cutting Length	Accepted Guide Wire	Handle	Unit	
PAP1-A2-20-35-OL 501028	Precut	2 Lumen, braided wire	20 mm	.035"	OL power connection	2	
PAP1-C2-20-35-OL 501034	Tapered	2 Lumen, braided wire	20 mm	.035"	OL power connection	2	
PAP1-C2-25-35-OL 501036	Tapered	2 Lumen, braided wire	25 mm	.035"	OL power connection	2	
PAP1-C2-30-35-OL 501039	Tapered	2 Lumen, braided wire	30 mm	.035"	OL power connection	2	
PAP1-C5-20-35-OL 501046	Tapered	2 Lumen, monofil	20 mm	.035"	OL power connection	2	
PAP1-C5-30-35-OL 501047	Tapered	2 Lumen, monofil	30 mm	.035"	OL power connection	2	
PAP1-E2-25-35-OL 501055	Long nose	2 Lumen, braided wire	25 mm	.035"	OL power connection	2	
PAP1-E2-30-35-OL 501056	Long nose	2 Lumen, braided wire	30 mm	.035"	OL power connection	2	
PAP1-K2-06-35-OL 501060	Needle	2 Lumen	6 mm	.035"	OL power connection	2	
PAP1-M2-20-35-OL 501840	B II (Type Erlangen)	2 Lumen, braided wire	20 mm	.035"	OL power connection	2	

= Marked products may have a delivery time of up to 8 weeks.

Errors excepted.

PAP1 _ Disposable papillotomes _ 3 Lumen




Product Code <i>Order Number</i>	Shape	Characteristics	Cutting Length	Accepted Guide Wire	Handle	Unit
PAP1-C7-20-35-OL 501048	Tapered	3 Lumen, braided wire	20 mm	.035"	OL power connection	2
PAP1-C7-30-35-OL 501050	Tapered	3 Lumen, braided wire	30 mm	.035"	OL power connection	2
PAP1-C8-20-35-OL 501051	Tapered	3 Lumen, monofil	20 mm	.035"	OL power connection	2
PAP1-C8-30-35-OL 501052	Tapered	3 Lumen, monofil	30 mm	.035"	OL power connection	2
PAP1-F8-30-35-OL 501993	Tapered	3 Lumen, monofil, insulated cutting wire (ProCut)	15 / 15 mm proximal insulated	.035"	OL power connection	2



PAP1 _ Disposable papillotomes for balloon enteroscopy (length 260 cm)








Product Code <i>Order Number</i>	Shape	Characteristics	Cutting Length	Accepted Guide Wire	Handle	Unit
PAP1-D2-30-21-OL 501853	Tapered	2 Lumen, braided wire	30 mm	.021"	OL power connection	2
PAP1-N2-20-21-OL 501902	B II (Type Erlangen)	2 Lumen, braided wire	20 mm	.021"	OL power connection	2



PAP Series Papillotomes

PAP1 _ Disposable **easycure** papillotomes



Product Code <i>Order Number</i>	Shape	Characteristics	Cutting Length	Accepted Guide Wire	Handle	Unit	
PAP1-A2-20-E35-OL 501994	Precut	3 Lumen, braided wire	20 mm	.035"	OL power connection	2	
PAP1-C7-30-E35-OL 501929	Tapered	3 Lumen, braided wire	30 mm	.035"	OL power connection	2	
PAP1-C8-30-E35-OL 501919	Tapered	3 Lumen, monofil	30 mm	.035"	OL power connection	2	
PAP1-F8-30-E35-OL 502017	Tapered	3 Lumen, monofil, insulated cutting wire (ProCut)	15 / 15 mm proximal insulated	.035"	OL power connection	2	
PAP1-K8-06-E35-OL 501920	Needle	3 Lumen	6 mm	.035"	OL power connection	2	

Errors excepted.