

medw^ork®


POL

Series

Polypectomy Snares


POL

POL^{Series} Polypectomy Snares



Polypectomy Snares

Numerous innovative models and the use of special cutting wires make this product really special. One hundred percent hand-made and „Made in Germany“ – top quality to the last detail.



Resection Master: POL1-C5-30-23-220-OL®

Benefit from the features:


- ideal for the lift-and-cut technique used for mucosa resection
- extremely high expansion and application force with monofilament flat wire
- loop has a rectangular design to prevent slippage

POL Series Polypectomy Snares

POL1 _ Disposable polypectomy snares _ Oval



Product Code <i>Order Number</i>	Shape	Characteristics	∅ Snare	∅ Tube	Length	Handle	Unit
POL1-B1-10-23-220-OL 501891	Oval	Braided wire	10 mm	2,3 mm	220 cm	OL power connection	10
POL1-B1-15-23-220-OL 501207	Oval	Braided wire	15 mm	2,3 mm	220 cm	OL power connection	10
POL1-B1-30-23-220-OL 501208	Oval	Braided wire	30 mm	2,3 mm	220 cm	OL power connection	10
POL1-B1-40-23-220-OL 501210	Oval	Braided wire	40 mm	2,3 mm	220 cm	OL power connection	10
POL1-B9-15-23-220-OL 502050	Oval	Monofil	15 mm	2,3 mm	220 cm	OL power connection	10
POL1-B9-22-23-220-OL 502048	Oval	Monofil	22 mm	2,3 mm	220 cm	OL power connection	10
POL1-B9-30-23-220-OL 502049	Oval	Monofil	30 mm	2,3 mm	220 cm	OL power connection	10
POL1-D1-30-23-220-OL 501223	Oval, combi	Braided wire	15 / 30 mm	2,3 mm	220 cm	OL power connection	10



POL1 _ Disposable polypectomy snares _ Oval, rotatable



Product Code <i>Order Number</i>	Shape	Characteristics	∅ Snare	∅ Tube	Length	Handle	Unit
POL1-B3-10-23-220-OL 502029	Oval	Monofil, rotatable	10 mm	2,3 mm	220 cm	OL power connection	10
POL1-B3-15-23-220-OL 501215	Oval	Monofil, rotatable	15 mm	2,3 mm	220 cm	OL power connection	10
POL1-B3-22-23-220-OL 501216	Oval	Monofil, rotatable	22 mm	2,3 mm	220 cm	OL power connection	10
POL1-B3-30-23-220-OL 501217	Oval	Monofil, rotatable	30 mm	2,3 mm	220 cm	OL power connection	10
POL1-B8-22-23-220-OL 501219	Oval	Braided wire, rotatable	22 mm	2,3 mm	220 cm	OL power connection	5
POL1-B8-30-23-220-OL 501220	Oval	Braided wire, rotatable	30 mm	2,3 mm	220 cm	OL power connection	5




POL Series Polypectomy Snares

POL1 _ Disposable polypectomy snares _ Oval, economy




Product Code Order Number	Shape	Characteristics	ø Snare	ø Tube	Length	Handle	Unit
POL1-B7-20-23-220-OL 501903	Oval	Braided wire	20 mm	2,3 mm	220 cm	OL power connection	10
POL1-B7-30-23-220-OL 501218	Oval	Braided wire	30 mm	2,3 mm	220 cm	OL power connection	10



POL1 _ Disposable polypectomy snares _ Combi instrument („Snare Inflator“)



Product Code Order Number	Shape	Characteristics	ø Snare	ø Tube	Length	Handle	Unit
POL1-H3-22-27-230-OL 501231	Combi instrument	Polypectomy Snare and injection needle in one instrument, monofil	22 mm	2,7 mm	230 cm	OL power connection	5
POL1-H3-30-27-230-OL 501948	Combi instrument	Polypectomy Snare and injection needle in one instrument, monofil	30 mm	2,7 mm	230 cm	OL power connection	5



POL1 _ Disposable polypectomy snares for mucosal resection („Resection Master“)




Product Code Order Number	Shape	Characteristics	ø Snare	ø Tube	Length	Handle	Unit
POL1-C5-22-23-220-OL 501896	Oval	Flat metal band, proximal strengthened	22 mm	2,3 mm	220 cm	OL power connection	5
POL1-C5-30-23-220-OL 501222	Oval	Flat metal band, proximal strengthened	30 mm	2,3 mm	220 cm	OL power connection	5



POL1 _ Disposable polypectomy snares_ Hexagonal



Product Code Order Number	Shape	Characteristics	ø Snare	ø Tube	Length	Handle	Unit
POL1-F1-30-23-220-OL 502101	Hexagonal	Braided wire	30 mm	2,3 mm	220 cm	OL power connection	10




POL Series Polypectomy Snares

POL1 _ Disposable polypectomy snares _ Asymmetrical, combi







Product Code Order Number	Shape	Characteristics	∅ Snare	∅ Tube	Length	Handle	Unit
POL1-K1-25-23-220-OL 502102	Asymmetrical, combi	Braided wire	15 / 25 mm	2,3 mm	220 cm	OL power connection	10



POL1 _ Disposable polypectomy snares for balloon enteroscopy




Product Code Order Number	Shape	Characteristics	∅ Snare	∅ Tube	Length	Handle	Unit
POL1-B1-30-18-260-OL 501843	Oval	Braided wire	30 mm	1,8 mm	260 cm	OL power connection	10
POL1-B1-30-23-260-OL 501209	 Oval	Braided wire	30 mm	2,3 mm	260 cm	OL power connection	10
POL1-E1-30-18-260-OL 501844	Oval, combi (Type Erlangen)	Braided wire, proximal strengthened by a flat spring	15 / 30 mm	1,8 mm	260 cm	OL power connection	10
POL1-E1-30-23-260-OL 501228	 Oval, combi (Type Erlangen)	Braided wire, proximal strengthened by a flat spring	15 / 30 mm	2,3 mm	260 cm	OL power connection	10

POL2 _ Reusable polypectomy snares _ Asymmetrical

134°C

Product Code Order Number	Shape	Characteristics	∅ Snare	∅ Tube	Length	Handle	Unit
POL2-A1-15-23-220 501320	Asymmetrical	Braided wire	15 mm	2,3 mm	220 cm	-	1
POL2-A1-25-23-220 501321	Asymmetrical	Braided wire	25 mm	2,3 mm	220 cm	-	1
POL2-A1-25-23-220-OL 501322	Asymmetrical	Braided wire	25 mm	2,3 mm	220 cm	OL power connection	1






 = Marked products may have a delivery time of up to 8 weeks.



Errors excepted.

POL Series Polypectomy Snares

POL2 _ Reusable polypectomy snares _ Oval

134°C

Product Code <i>Order Number</i>	Shape	Characteristics	∅ Snare	∅ Tube	Length	Handle	Unit
POL2-B1-15-23-220 501338	 Oval	Braided wire	15 mm	2,3 mm	220 cm	-	1
POL2-B1-30-23-220 501339	Oval	Braided wire	30 mm	2,3 mm	220 cm	-	1
POL2-B1-30-23-220-OL 501340	Oval	Braided wire	30 mm	2,3 mm	220 cm	OL power connection	1
POL2-B1-40-23-220 501342	Oval	Braided wire	40 mm	2,3 mm	220 cm	-	1
POL2-B1-40-23-220-OL 501343	Oval	Braided wire	40 mm	2,3 mm	220 cm	OL power connection	1
POL2-B2-30-23-220-OL 501351	Oval	Braided wire, insulated tip	30 mm	2,3 mm	220 cm	OL power connection	1
POL2-B2-40-23-220-OL 501353	Oval	Braided wire, insulated tip	40 mm	2,3 mm	220 cm	OL power connection	1
POL2-B9-15-23-220-OL 502113	Oval	Monofil	15 mm	2,3 mm	220 cm	OL power connection	1
POL2-B9-20-23-220 502116	 Oval	Monofil	20 mm	2,3 mm	220 cm	-	1
POL2-B9-30-23-220 502115	 Oval	Monofil	30 mm	2,3 mm	220 cm	-	1
POL2-B9-30-23-220-OL 502114	Oval	Monofil	30 mm	2,3 mm	220 cm	OL power connection	1



 = Marked products may have a delivery time of up to 8 weeks.


Errors excepted.

POL Series Polypectomy Snares

POL2 _ Reusable polypectomy snares _ Oval, rotatable

134°C


Product Code Order Number	Shape	Characteristics	∅ Snare	∅ Tube	Length	Handle	Unit
POL2-B8-22-23-220-OL 501845	Oval	Braided wire, rotatable	22 mm	2,3 mm	220 cm	OL power connection	1
POL2-B8-30-23-220-OL 501366	Oval	Braided wire, rotatable	30 mm	2,3 mm	220 cm	OL power connection	1



POL2 _ Reusable polypectomy snares _ Oval (Type Erlangen)

134°C

Product Code Order Number	Shape	Characteristics	∅ Snare	∅ Tube	Length	Handle	Unit
POL2-C1-15-23-220 502147	 Oval (Type Erlangen)	Braided wire, proximal strengthened by a flat spring	15 mm	2,3 mm	220 cm	-	1
POL2-C1-20-23-220 502148	 Oval (Type Erlangen)	Braided wire, proximal strengthened by a flat spring	20 mm	2,3 mm	220 cm	-	1
POL2-C1-30-23-220 501368	Oval (Type Erlangen)	Braided wire, proximal strengthened by a flat spring	30 mm	2,3 mm	220 cm	-	1
POL2-C1-30-23-220-OL 501369	Oval (Type Erlangen)	Braided wire, proximal strengthened by a flat spring	30 mm	2,3 mm	220 cm	OL power connection	1
POL2-C1-50-23-220 501370	Oval (Type Erlangen)	Braided wire, proximal strengthened by a flat spring	50 mm	2,3 mm	220 cm	-	1
POL2-C1-50-23-220-OL 501371	Oval (Type Erlangen)	Braided wire, proximal strengthened by a flat spring	50 mm	2,3 mm	220 cm	OL power connection	1
POL2-C2-30-23-220 501375	Oval (Type Erlangen)	Braided wire, proximal strengthened by a flat spring, insulated tip	30 mm	2,3 mm	220 cm	-	1
POL2-C2-30-23-220-OL 501376	Oval (Type Erlangen)	Braided wire, proximal strengthened by a flat spring, insulated tip	30 mm	2,3 mm	220 cm	OL power connection	1
POL2-C2-50-23-220 501378	Oval (Type Erlangen)	Braided wire, proximal strengthened by a flat spring, insulated tip	50 mm	2,3 mm	220 cm	-	1
POL2-C2-50-23-220-OL 501379	Oval (Type Erlangen)	Braided wire, proximal strengthened by a flat spring, insulated tip	50 mm	2,3 mm	220 cm	OL power connection	1



 = Marked products may have a delivery time of up to 8 weeks.


Errors excepted.

POL Series Polypectomy Snares

POL2 _ Reusable polypectomy snares _ Oval, combi (type Erlangen)

134°C


Product Code Order Number	Shape	Characteristics	∅ Snare	∅ Tube	Length	Handle	Unit
POL2-E1-30-23-220 501404	Oval, combi (Type Erlangen)	Braided wire, proximal strengthened by a flat spring	15 / 30 mm	2,3 mm	220 cm	-	1
POL2-E1-30-23-220-OL 501405	Oval, combi (Type Erlangen)	Braided wire, proximal strengthened by a flat spring	15 / 30 mm	2,3 mm	220 cm	OL power connection	1



POL2 _ Reusable polypectomy snares_ Hexagonal

134°C

Product Code Order Number	Shape	Characteristics	∅ Snare	∅ Tube	Length	Handle	Unit
POL2-F1-30-23-220 501409	Hexagonal	Braided wire	30 mm	2,3 mm	220 cm	-	1
POL2-F1-30-23-220-OL 501410	Hexagonal	Braided wire	30 mm	2,3 mm	220 cm	OL power connection	1



HAN2 _ Reusable accessories for the POL series

134°C

Product Code Order Number	Characteristics	Lift	Unit
HAN2-100-OL 500652	Reusable handle, rotatable, OL power connection	100 mm	1
HAN2-100-ST 500653	Reusable handle, rotatable, standard power connection	100 mm	1
HAN2-Adapter-OL-ST 500654	Reusable adapter, converts OL into standard power connection	-	1
MIS2-HF-05-OL 500962	Reusable high-frequency-cable, 5 mm plug on the part of appliance, OL HF power connector	-	1

