


medw<sup>o</sup>rk®

# PRO

Series

**Plastic Stents**

PRO



## Plastic ERCP stents

We offer a wide range of products that all have one thing in common: the optimum combination of stiffness and flexibility, ideal for the most demanding of tasks. You have the choice of conventional or completely preloaded insertion kits for the safe and efficient placement of plastic stents.



### Complete sets: PRO3-A3-10-10-12<sup>®</sup>

The complete set comprises of the guide catheter, pusher and plastic stent makes effective and atraumatic placement of the stent the bile duct a simple task.

#### Benefit from the features:

- \_ ready for immediate use without unnecessary loss of time
- \_ makes your storage management much simpler
- \_ colour-coded system components

# PRO Series Plastic Stents

## PRO1 \_ Disposable Endoprosthesis \_ 7 F, straight, 2 flaps, bile duct, hard, tapered



Product Code <i>Order Number</i>	Length between Flaps / Total	Unit
PRO1-A1-07-04-06 501461	4 cm / 6 cm	1
PRO1-A1-07-06-08 501463	6 cm / 8 cm	1

Product Code <i>Order Number</i>	Length between Flaps / Total	Unit
PRO1-A1-07-09-11 501466	9 cm / 11 cm	1
PRO1-A1-07-12-14 501469	12 cm / 14 cm	1



Intermediate lengths available on request.

## PRO1 \_ Disposable Endoprosthesis \_ 8.5 F, straight, 2 flaps, bile duct, hard, tapered



Product Code <i>Order Number</i>	Length between Flaps / Total	Unit
PRO1-A1-08-04-06 501474	4 cm / 6 cm	1
PRO1-A1-08-06-08 501476	6 cm / 8 cm	1

Product Code <i>Order Number</i>	Length between Flaps / Total	Unit
PRO1-A1-08-09-11 501479	9 cm / 11 cm	1
PRO1-A1-08-12-14 501482	12 cm / 14 cm	1



Intermediate lengths available on request.

## PRO1 \_ Disposable Endoprosthesis \_ 10 F, straight, 2 flaps, bile duct, hard, tapered



Product Code <i>Order Number</i>	Length between Flaps / Total	Unit
PRO1-A1-10-04-06 501489	4 cm / 6 cm	1
PRO1-A1-10-06-08 501491	6 cm / 8 cm	1

Product Code <i>Order Number</i>	Length between Flaps / Total	Unit
PRO1-A1-10-09-11 501494	9 cm / 11 cm	1
PRO1-A1-10-12-14 501497	12 cm / 14 cm	1



Intermediate lengths available on request.

## PRO1 \_ Disposable Endoprosthesis \_ 11.5 F, straight, 2 flaps, bile duct, hard, tapered



Product Code <i>Order Number</i>	Length between Flaps / Total	Unit
PRO1-A1-11-04-06 501504	4 cm / 6 cm	1
PRO1-A1-11-06-08 501506	6 cm / 8 cm	1

Product Code <i>Order Number</i>	Length between Flaps / Total	Unit
PRO1-A1-11-09-11 501509	9 cm / 11 cm	1
PRO1-A1-11-12-14 501512	12 cm / 14 cm	1



Intermediate lengths available on request.

Errors excepted.

## PRO1 \_ Disposable Endoprosthesis \_ 8.5 F, straight, 8 flaps, bile duct, hard, tapered



Product Code <i>Order Number</i>	Length between Flaps / Total	Unit
PRO1-A2-08-04-06 501515	4 cm / 6 cm	1
PRO1-A2-08-06-08 501516	6 cm / 8 cm	1

Product Code <i>Order Number</i>	Length between Flaps / Total	Unit
PRO1-A2-08-09-11 501517	9 cm / 11 cm	1
PRO1-A2-08-12-14 501518	12 cm / 14 cm	1



Intermediate lengths available on request.

## PRO1 \_ Disposable Endoprosthesis \_ 10 F, straight, 8 flaps, bile duct, hard, tapered



Product Code <i>Order Number</i>	Length between Flaps / Total	Unit
PRO1-A2-10-04-06 501519	4 cm / 6 cm	1
PRO1-A2-10-06-08 501521	6 cm / 8 cm	1

Product Code <i>Order Number</i>	Length between Flaps / Total	Unit
PRO1-A2-10-09-11 501524	9 cm / 11 cm	1
PRO1-A2-10-12-14 501527	12 cm / 14 cm	1



Intermediate lengths available on request.

## PRO1 \_ Disposable Endoprosthesis \_ 11.5 F, straight, 8 flaps, bile duct, hard, tapered



Product Code <i>Order Number</i>	Length between Flaps / Total	Unit
PRO1-A2-11-04-06 501531	4 cm / 6 cm	1
PRO1-A2-11-06-08 501532	6 cm / 8 cm	1

Product Code <i>Order Number</i>	Length between Flaps / Total	Unit
PRO1-A2-11-09-11 501534	9 cm / 11 cm	1
PRO1-A2-11-12-14 501535	12 cm / 14 cm	1



Intermediate lengths available on request.

# PRO Series Plastic Stents

## PRO1 \_ Disposable Endoprosthesis \_ 7 F, straight, 2 flaps, bile duct, soft, tapered



Product Code <i>Order Number</i>	Length between Flaps / Total	Unit
PRO1-A3-07-04-06 501859	4 cm / 6 cm	1
PRO1-A3-07-06-08 501863	6 cm / 8 cm	1

Product Code <i>Order Number</i>	Length between Flaps / Total	Unit
PRO1-A3-07-09-11 501867	9 cm / 11 cm	1
PRO1-A3-07-12-14 501871	12 cm / 14 cm	1



Intermediate lengths available on request.

## PRO1 \_ Disposable Endoprosthesis \_ 8.5 F, straight, 2 flaps, bile duct, soft, tapered



Product Code <i>Order Number</i>	Length between Flaps / Total	Unit
PRO1-A3-08-04-06 501860	4 cm / 6 cm	1
PRO1-A3-08-06-08 501864	6 cm / 8 cm	1

Product Code <i>Order Number</i>	Length between Flaps / Total	Unit
PRO1-A3-08-09-11 501868	9 cm / 11 cm	1
PRO1-A3-08-12-14 501872	12 cm / 14 cm	1



Intermediate lengths available on request.

## PRO1 \_ Disposable Endoprosthesis \_ 10 F, straight, 2 flaps, bile duct, soft, tapered



Product Code <i>Order Number</i>	Length between Flaps / Total	Unit
PRO1-A3-10-04-06 501861	4 cm / 6 cm	1
PRO1-A3-10-06-08 501865	6 cm / 8 cm	1

Product Code <i>Order Number</i>	Length between Flaps / Total	Unit
PRO1-A3-10-09-11 501869	9 cm / 11 cm	1
PRO1-A3-10-12-14 501873	12 cm / 14 cm	1



Intermediate lengths available on request.

## PRO1 \_ Disposable Endoprosthesis \_ 11.5 F, straight, 2 flaps, bile duct, soft, tapered



Product Code <i>Order Number</i>	Length between Flaps / Total	Unit
PRO1-A3-11-04-06 501862	4 cm / 6 cm	1
PRO1-A3-11-06-08 501866	6 cm / 8 cm	1

Product Code <i>Order Number</i>	Length between Flaps / Total	Unit
PRO1-A3-11-09-11 501870	9 cm / 11 cm	1
PRO1-A3-11-12-14 501874	12 cm / 14 cm	1



Intermediate lengths available on request.

Errors excepted.

## PRO1 \_ Disposable endoprosthesis \_ Pigtail, double, bile duct




Product Code <i>Order Number</i>	Shape	Characteristics	Diameter	Length between Flaps	Total Length	Unit
PR01-D4-08-06-09 501580	Pigtail, double	Punched	8,5 F	6 cm	9 cm	1
PR01-D4-08-09-12 501581	Pigtail, double	Punched	8,5 F	9 cm	12 cm	1
PR01-D4-08-12-15 501582	Pigtail, double	Punched	8,5 F	12 cm	15 cm	1
PR01-D4-10-06-09 501583	Pigtail, double	Punched	10 F	6 cm	9 cm	1
PR01-D4-10-09-12 501584	Pigtail, double	Punched	10 F	9 cm	12 cm	1
PR01-D4-10-12-15 501585	Pigtail, double	Punched	10 F	12 cm	15 cm	1



## PRO1 \_ Disposable endoprosthesis insertion sets



Product Code <i>Order Number</i>	Characteristics	∅ Guiding Catheter	Length	Accepted Guide Wire	Unit
PR01-H7-05-205-32 502124	Radiopaque tip, LL connection, for 8.5 F endoprosthesis	5 F	205 cm	.032"	1
PR01-H7-06-205-35 501908	Radiopaque tip, LL connection, for 10 F endoprosthesis	6 F	205 cm	.035"	1
PR01-H7-07-205-35 501909	Radiopaque tip, LL connection, for 11.5 F endoprosthesis	7 F	205 cm	.035"	1



## PRO1 \_ Disposable **easycure** endoprosthesis insertion sets



Product Code <i>Order Number</i>	Characteristics	∅ Guiding Catheter	Length	Accepted Guide Wire	Unit
PR01-H7-05-205-E32 502123	Radiopaque tip, LL connection, for 8.5 F endoprosthesis	5 F	205 cm	.032"	1
PR01-H7-06-205-E35 501918	Radiopaque tip, LL connection, for 10 F endoprosthesis	6 F	205 cm	.035"	1



# PRO Series Plastic Stents

## PRO1 \_ Disposable guiding catheters



Product Code <i>Order Number</i>	Characteristics	Diameter	Length	Accepted Guide Wire	Unit
PR01-F5-05-220-32 502119	For 8.5 F endoprosthesis, distal radiopaque marking, white	5 F	220 cm	.032"	1
PR01-F5-05-320-32 502118	For 8.5 F endoprosthesis, distal radiopaque marking, white	5 F	320 cm	.032"	1
PR01-F5-06-220-35 502105	For 10 F endoprosthesis, distal radiopaque marking, blue	6 F	220 cm	.035"	1
PR01-F5-06-320-35 502106	For 10 F endoprosthesis, distal radiopaque marking, blue	6 F	320 cm	.035"	1
PR01-F5-07-220-35 502107	For 11.5 F endoprosthesis, distal radiopaque marking, white	7 F	220 cm	.035"	1
PR01-F5-07-320-35 502108	For 11.5 F endoprosthesis, distal radiopaque marking, white	7 F	320 cm	.035"	1

## PRO1 \_ Disposable pusher




Product Code <i>Order Number</i>	Characteristics	Diameter	Length	Unit
PR01-G7-07-170 502109	Radiopaque	7 F	170 cm	1
PR01-G7-08-170 502110	Radiopaque	8,5 F	170 cm	1
PR01-G7-10-170 502111	Radiopaque	10 F	170 cm	1
PR01-G7-11-170 502112	Radiopaque	11,5 F	170 cm	1

## PRO3 \_ Disposable endoprosthesis insertion sets with pre-loaded 8.5 F plastic stent (straight, 2 flaps, soft, for .032" guiding catheter)



Product Code Order Number	Length between Flaps / Total	Unit
PR03-A3-08-04-06 502083	4 cm / 6 cm	1
PR03-A3-08-05-07 502084	5 cm / 7 cm	1
PR03-A3-08-06-08 502085	6 cm / 8 cm	1
PR03-A3-08-07-09 502086	7 cm / 9 cm	1
PR03-A3-08-08-10 502087	8 cm / 10 cm	1

Product Code Order Number	Length between Flaps / Total	Unit
PR03-A3-08-09-11 502088	9 cm / 11 cm	1
PR03-A3-08-10-12 502089	10 cm / 12 cm	1
PR03-A3-08-11-13 502090	11 cm / 13 cm	1
PR03-A3-08-12-14 502091	12 cm / 14 cm	1




## PRO3 \_ Disposable endoprosthesis insertion sets with pre-loaded 10 F plastic stent (straight, 2 flaps, soft, for .032" guiding catheter)



Product Code Order Number	Length between Flaps / Total	Unit
PR03-A3-10-04-06 502061	4 cm / 6 cm	1
PR03-A3-10-05-07 502060	5 cm / 7 cm	1
PR03-A3-10-06-08 502053	6 cm / 8 cm	1
PR03-A3-10-07-09 502059	7 cm / 9 cm	1
PR03-A3-10-08-10 502058	8 cm / 10 cm	1

Product Code Order Number	Length between Flaps / Total	Unit
PR03-A3-10-09-11 502057	9 cm / 11 cm	1
PR03-A3-10-10-12 502056	10 cm / 12 cm	1
PR03-A3-10-11-13 502055	11 cm / 13 cm	1
PR03-A3-10-12-14 502054	12 cm / 14 cm	1



## PRO3 \_ Disposable endoprosthesis insertion sets with pre-loaded 11.5 F plastic stent (straight, 2 flaps, soft, for .032" guiding catheter)



Product Code Order Number	Length between Flaps / Total	Unit
PR03-A3-11-04-06 502092	4 cm / 6 cm	1
PR03-A3-11-05-07 502093	5 cm / 7 cm	1
PR03-A3-11-06-08 502094	6 cm / 8 cm	1
PR03-A3-11-07-09 502095	7 cm / 9 cm	1
PR03-A3-11-08-10 502096	8 cm / 10 cm	1

Product Code Order Number	Length between Flaps / Total	Unit
PR03-A3-11-09-11 502097	9 cm / 11 cm	1
PR03-A3-11-10-12 502098	10 cm / 12 cm	1
PR03-A3-11-11-13 502099	11 cm / 13 cm	1
PR03-A3-11-12-14 502100	12 cm / 14 cm	1

