


medw^ork®

WIR

Series

Guide Wires



Guide Wires

Copes with even difficult stenosis and provides atraumatic access to the bile duct with products of the WIR Series.

gSlider: WIR1-F9-35-450[®]

Benefit from the features:

- _ hydrophilic tip provides revolutionary sliding properties with the innovative Aquarius™ coating
- _ excellent radiological visibility
- _ maximum guide wire control
- _ extremely kink resistant tip

WIR^{Series} Guide Wires

WIR1 _ Disposable guide wires _ Nitinol, hydrophilic tip („gSlider“)



Product Code <i>Order Number</i>	Shape	Characteristics	Diameter	Length	Unit	
WIR1-F9-35-260 502145	Straight	Nitinol, insulated, marked, green/white, radiopaque hydrophilic tip	.035"	260 cm	2	
WIR1-F9-35-450 502144	Straight	Nitinol, insulated, marked, green/white, radiopaque hydrophilic tip	.035"	450 cm	2	

WIR1 _ Disposable guide wires _ Steel



Product Code <i>Order Number</i>	Shape	Characteristics	Diameter	Length	Unit	
WIR1-A2-35-110 501773	Straight	Steel, coated	.035"	110 cm	2	
WIR1-A2-35-400 501775	Straight	Steel, coated	.035"	400 cm	2	

WIR1 _ Disposable guide wires _ Nitinol



Product Code <i>Order Number</i>	Shape	Characteristics	Diameter	Length	Unit	
WIR1-F1-18-400 501781	Straight	Nitinol, coated, radiopaque tip	.018"	400 cm	2	
WIR1-F1-21-400 501782	Straight	Nitinol, coated, radiopaque tip	.021"	400 cm	2	
WIR1-F1-21-500 501911	Straight	Nitinol, coated, radiopaque tip	.021"	500 cm	2	
WIR1-F5-25-400 501789	Straight	Nitinol, insulated, marked, black/white	.025"	400 cm	2	
WIR1-F5-35-400 501792	Straight	Nitinol, insulated, marked, black/white	.035"	400 cm	2	
WIR1-F7-35-260 501837	Straight	Nitinol, coated, extra stiff, radiopaque tip	.035"	260 cm	2	
WIR1-F7-35-400 501838	Straight	Nitinol, coated, extra stiff, radiopaque tip	.035"	400 cm	2	
WIR1-G5-35-400 501804	Straight	Nitinol, insulated, marked, black/white, extra soft	.035"	400 cm	2	

Errors excepted.

WIR1 _ Disposable guide wires _ Eder-Püstow



Product Code <i>Order Number</i>	Shape	Characteristics	Diameter	Length	Unit	
WIR1-E1-72-260 501780	Straight	Steel, flexible tip, Eder-Püstow	.072"	260 cm	2	
WIR1-E1-72-315 501783	Straight	Steel, flexible tip, Eder-Püstow	.072"	315 cm	2	

WIR1 _ Disposable guide wires _ Nitinol, J-tip



Product Code <i>Order Number</i>	Shape	Characteristics	Diameter	Length	Unit	
WIR1-H5-35-400 501806	J-Tip	Nitinol, insulated, marked, black/white	.035"	400 cm	2	

WIR1 _ Disposable **easyCure** guide wires



Product Code <i>Order Number</i>	Shape	Characteristics	Diameter	Length	Unit	
WIR1-F1-E21-200 501922	Straight	Nitinol, coated, radiopaque tip	.021"	200 cm	2	
WIR1-F5-E35-200 501921	Straight	Nitinol, insulated, marked, black/white	.035"	200 cm	2	

= Marked products may have a delivery time of up to 8 weeks.

Errors excepted.